
Instruction Manual

Instruction Manual.....	2
Manual de Instrucciones	12
Manuel d'instructions.....	23
Bedienungsanleitung	34

To immediately shut off power to Vortexer device, pull out the AC cord from the lower central back area of the device

Thank you for purchasing the Heathrow Scientific Vortexer. Your Vortexer mixing instrument is designed for general laboratory use

Safety Precautions

Your Vortexer instrument is designed for general laboratory use. This instruction manual cannot address all safety hazards. It is the responsibility of the user to consult and observe all health and safety precautions and to assess the instrument's suitability to perform the task.

Caution!

Before starting the Vortexer for the first time, please read this instruction manual carefully. This manual contains important information, safeguards, and operating instructions.

1. Never use this product in any manner inconsistent with these instructions.
2. This product is intended to be used indoors only.
3. Always ensure that the disconnecting plug (power plug) is accessible.
4. This product is intended for vortexing and mixing aqueous solutions in approved test tubes compatible to the speed and power of this unit. Note on Vortexing: use only approved containers with screw-tight or lockable lids.
5. The product is not intended to vortex very dense materials. Do not use materials with a density greater than 1.2g/cm.
6. Observe general laboratory safety precautions and regulations when using this product.
7. Do not attempt to operate with the mixing head attachment removed. Vortexer may become destabilized in operation.
8. The mixing head attachment must always be securely fastened when in operation: Always make certain that the head attachment is seated completely onto the top shaft before operating.
9. Do not use if the Vortexer is damaged.
10. Do not put hands onto or under the mixing head attachment unless the Vortexer is completely stopped.
11. Do not move unit while the the Vortexer is in operation.
12. Do not immerse the product in water. Injury due to shock or fire may result.
13. The mixing head attachment must be loaded symmetrically. Operating this unit with an unbalanced load will cause severe vibration and damage, which is not covered by warranty.
14. Do not fill tubes while inserted onto the Vortexer mixing head. Spilled liquid may harm the unit or the user.
15. Close tube lids/caps and confirm that they are completely secure before starting the Vortexer instrument. Open tube lids can be torn off during vortexing and damage the unit. Open lids can cause fluid to be dispersed inside the vortex instrument.
16. Do not insert non-standard tubes or other foreign articles into the Vortexer head.
17. Follow all precautions from the Material Safety Data Sheet (MSDS) for any reagent you use with this instrument.
18. Do not use solvent or flammable liquids near this or other electrical equipment.
19. Do not vortex materials that in combination or alone are flammable, explosive, or corrosive or may cause any other hazard.
20. Do not operate in a hazardous or flammable environment.

21. Do not use any solvent on the unit that may degrade the plastic of the Vortexer main housing or head.
22. Ensure the mixing head is protected from corrosion and mechanical damage. The mixing head must be cleaned with a pH-neutral cleaning liquid.
23. Always work in a manner which endangers neither the user nor any other person.
24. Should the instrument fail to work to specification, immediately stop vortexing. Contact the authorized dealer or the manufacturer if trouble persists.
25. Do not attempt to make contact with the mixing head attachment while unit is running. Doing so may cause the unit to fail and will void the warranty.
26. Only use original manufacturer's mixing head attachment, AC adapters, and any other spare parts.
27. If unit is malfunctioning or needs to be returned, please contact either the authorized dealer or Heathrow Scientific® LLC .
28. Improper use of the Vortexer instrument will void the warranty. If there is a failure during the warranty period, contact Heathrow Scientific® LLC for warranty service.
29. Do not attempt to use units that have not been correctly installed or repaired.
30. Do not attempt to disassemble or modify this product. Any attempt to do so will void the warranty.
31. Store at room temperature in a dry area. Do not expose to sunlight, moisture, or extreme temperatures for prolonged periods of time.
32. Before using any cleaning or decontamination methods please contact Heathrow Scientific® LLC to ensure that the proposed method will not damage the equipment. A pH-neutral cleaning liquid must be used.
33. Vortexer must be unplugged from electrical source before any cleaning or maintenance is performed.

To immediately shut off power to Vortexer device, pull out the AC cord from the lower central back area of the device

Attention: If main electrical power is interrupted, Vortexer device will immediately begin operating if 3-position switch has not been set to **OFF** position when power is restored.

If power is interrupted: 1) Unplug cord from back of device. 2) Set Vortexer switch to **OFF**

Warning!

- Always check for device movement when running in continuous mode. If movement of device is noted, adjust speed accordingly. Different loading of mixing head may result in reaching a resonance frequency whereby the device will move on the work surface/ table top.
- If the unit is not used as recommended by the manufacturer, the overall safety will be impaired.
- Check MSDS, wear required Personal Protective Equipment, and observe all applicable local and national regulations before dispensing and disposing of hazardous samples.
- The laboratory procedure should ensure that no person or hazardous substances are present within a 12" (30cm) zone around the Vortexer instrument when operating.
- This product does not contain bio-seals as per IEC/EN/CSA 61010-2-20 and cannot provide any level of containment in case of a spill or release of toxic, radioactive, or pathogenic micro-organisms thus these materials are not recommended to be used in this product.

Specifications

Product Description/Features: This device is a benchtop laboratory Vortexer with a variable speed control (OFF, and up to approx. **3000 rpm**), capable of vortexing microplates and 6 different tube sizes with the convenience of a single mixing head attachment. This Vortexer has enough flexibility to perform the standard tasks that a vortex mixer is required to perform, plus it has the ability to hold a variety of tubes, containers, and microplates all within one mixing head design.

Supplied Equipment and Spare Parts:

- 1 ea. Vortexer (classified as Installation Category 2)
 - HS120209, US/North America: 120VAC (110V – 120V, 50-60 HZ, 55 Watts, 1.0 Amps)
 - HS120318, Japan (PSE cord): 100VAC (100V – 110V, 50-60 Hz, 66 Watts, 1.3 Amps)
 - HS120214, United Kingdom: 230V (220V – 240V, 50-60 Hz, 44 Watts, .4 Amps)
 - HS120212, EU: 230VAC (220V – 240V, 50-60 Hz, 44 Watts, .4 Amps)
 - HS120210, Australian: 230VAC (220V – 240V, 50-60 Hz, 44 Watts, .4 Amps)
- 1 ea. AC power cord (WARNING only use power cords supplied by Heathrow Scientific):
 - HS120329: N. America
 - HS120488: Japan PSE rated
 - HS120331: UK
 - HS120332: CE
 - HS120333: Australia
- 1 ea. mixing head, HS120327
- 3 additional elastic retention cords, HS120422
- 3 additional central pads for intermittent mixing area, HS120328

Technical Features:

Speed range:

From 0 to approximately 3000 rpm - shaking to vortexing

Operation:

1) 3-position slider power switch:

Top position: Intermittent operation

Center position: OFF

Bottom position: Continuously ON

2) Infinite position rotary speed control for speeds from Off up to approximately 3000 rpm

Dimensions: 6.75" W x 7.5" D x 6" Height of main unit. 7.75" Height with Mixing Head attached. (171mm W x 190mm D x 197mm H)

Weight: 8.66 lbs (3.93 kg) – total Vortexer weight with Mixing Head attached

Orbital Diameter: 3.7mm

Electrical: See HS part numbers listed above.

Note: Product has a Thermal Protection Device.

Environmental:

For Indoor Use Only

Altitude: up to 2,000 m

Operating Temperature: 2° C to 40° C (35.6° F to 104° F), Max. Humidity 95%

Max. Relative Humidity ≤ 80% non-condensing

This mark is the confirmation that the unit conforms to the EU guidelines and has been tested according to the following EU Declaration of Conformity.

Conforms to:
IEC/EN 61010-1:2010, ed. 3.1, 2017

IEC/EN 61010-2-051, ed. 3, 2016

UL Std. 61010-1, ed. 3, 2016

710213

UL Std. 61010-2-051, ed. 3 (R2016)

CAN/CSA C22.2 No. 61010-1-12 (R2017)

CAN/CSA C22.6 No. 61010-2-051 (R2016)

This Product is in Compliance with: RoHS 3, PSE Certified Cord (Japan), WEEE

EC Declaration Of Conformity

Manufacturer: Heathrow Scientific® LLC, 620 Lakeview Parkway, Vernon Hills, IL 60061 USA

European Contact: Emergo, Europe, Princessegracht 20, 2514AP, The Hague, The Netherlands

Product: Vortexer

Product Model No(s): 120209, 120210, 120212, 120214, 120318

This unit has been constructed and conforms to the following:

EMC Standards and Directives:

EN 61326-1:2013 (IEC 61326-1:2012)	2014/35/EU (Low Voltage Directive)
FCC Part 15, Subpart B: 2015	2012/19/EU (WEEE)
ICES-003, Issue 6	2015/863/EU (RoHS3)

Safety Standards:

IEC/EN 61010-1:2010, ed. 3.1, 2017
IEC/EN 61010-2-051, ed. 3, 2016
UL Std. 61010-1, ed. 3, 2016
UL Std. 61010-2-051, ed. 3 (R2016)
CAN/CSA C22.2 No. 61010-1-12 (R2017)
CAN/CSA C22.6 No. 61010-2-051 (R2016)

Signed:

Gary Kamees

Title: Vice President Product Development/Manufacturing

Location: Vernon Hills, IL

Date of Issue: June 2022

UKCA Declaration Of Conformity

Manufacturer: Heathrow Scientific® LLC, 620 Lakeview Parkway, Vernon Hills, IL 60061 USA

Product: Vortexer

Product Model No(s): 120209, 120210, 120212, 120214, 120318

This unit has been constructed and conforms to the following:

Directives and Regulations:

Electrical Equipment (Safety) Regulations 2016
The Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012
The Electromagnetic Compatibility Regulations 2016

EMC Standards:

EN 61326-1:2013 (IEC61326-1:2012)
ICES-003, Issue 6

Safety Standards:

IEC/EN 61010-1:2010, ed. 3.1, 2017	IEC/EN 61010-2-051, ed. 3, 2016
UL Std. 61010-1, ed. 3, 2016	UL Std. 61010-2-051, ed. 3 (R2016)
CAN/CSA C22.2 No. 61010-1-12 (R2017)	CAN/CSA C22.6 No. 61010-2-051 (R2016)

Signed:

Gary Kamees

Title: Vice President Product Development/Manufacturing


Location: Vernon Hills, IL

Date of Issue: October 2022

Your Vortexer is supplied with a pop-on mixing head attachment and an AC power cord located in the top corrugated box insert **1** of the product packaging. After you fully unpack your Vortexer from the package and remove the protective plastic bag, install the mixing head attachment (check that it is right side up).


IMPORTANT: YOU MUST PULL OUT INNER FRONT SUPPORT CARD **3** BEFORE REMOVING YOUR HEATHROW SCIENTIFIC® VORTEXER. USE TWO HANDS TO GRASP THE UNIT WHEN LIFTING OUT OF THE BOX. **4**

TO INSTALL THE MIXING HEAD ATTACHMENT YOU MUST PRESS DOWN FIRMLY WITH BOTH HANDS UNTIL YOU FEEL THE MIXING HEAD ATTACHMENT SNAP FULLY INTO PLACE.

This will ensure that the mixing head attachment will not come off the Vortexer during operation.

Place the Vortex Mixer on a level, stable surface near the grounded electrical outlet. The surface should be clean and free of dust to ensure that the feet grip the surface firmly. Allow sufficient clearance on all sides of the unit for proper ventilation and usage. With the power switch in the **OFF** position (central on the vertical 3-position), plug the power cord into the back of the Vortex Mixer.


Apply even pressure to ensure the head is secure

Note: On units equipped with a frequency switch (located above the power cord connection on rear) set switch to either 50Hz or 60Hz to match the local power company supply. Turn the speed rotary dial turned all the way to the left- (in the OFF position). Plug the device into a properly grounded outlet.

Power switch in the **OFF** position - center

50Hz-60Hz frequency switch location

Back view of Vortexer

To immediately shut off power to Vortexer device, pull out the AC cord from the lower central back area of the device

Modes of Operation:


WARNING: If the Vortexer is used in a manner not specified by Heathrow Scientific, the protection provided by the equipment may be impaired.

Duty Cycle for ALL VERSIONS: Maximum run time of 20 minutes followed by 60 minute OFF time when in Continuous Operation mode.

100V-110V Vortexer: Set power input frequency switch on rear of unit to either 50Hz or 60Hz to match local power company supply.

220V-240V Vortexer: Set power input frequency switch on rear of unit to either 50Hz or 60Hz to match local power company supply.


The rotary speed control dial, located in the lower center of the front panel, has infinite adjustability (**OFF**, and **0 to approx 3000 rpm**). Vortexing speed increases to a maximum of approx. 3,000rpm as the knob is turned clockwise. Before turning on the Vortexer, set the speed dial to **OFF** position (rotate counter-clockwise until the knob stops and a “click” is felt). To begin operation, move the power switch to either the Intermittent **TOUCH** or **ON** position and turn the speed control to the desired setting. If the power switch is in the **ON** position, the mixing head will begin moving immediately. If the power switch is in the Intermittent **TOUCH** position, operation will begin when the rotary speed control dial is rotated clockwise to the maximum speed and an object (test tube) is placed in the center of the mixing head, held at the top of the test tube, and moderate pressure is applied to push down the mixing head. Mixing will cease when the object is removed and pressure is released.

I - Intermittent - Top position
II - Continuous - Bottom position


Intermittent/touch mode:

Place the 3-position switch in the Top Intermittent **TOUCH** position and rotate the speed control dial clockwise to maximum speed. When pressure is released from the mixing head, operation will stop. Intermittent operation will begin when a test tube is pressed down in the center area of the mixing head, while the tube is held firmly at the top.


Hands free/continuous on mode:

Place the 3-position switch in the **ON** (bottom) position and set the rotary speed control as described above. Hands Free mode allows mixing or vortexing of the tubes or microplate and will continue to run until you stop the mixer by placing the 3-position switch to the **OFF** position.

Note: To prevent overheating of the Vortexer during continuous use, there is an internal Temperature Protection Device that will temporarily shut off the device. Once the device cools down, it will be operational again. Please set switch to **OFF** position whenever it is not in operation.

Attention: It is recommended to **NOT** leave the Vortexer unattended while in operation.

OFF mode:


Always place the 3-position switch to the **OFF** position (center) when not in use!


Loading Tubes

The Vortexer mixing head is designed to hold several different tube sizes for hands-free mixing/vortexing.

Caution: Always load onto mixing head so that the tubes are balanced



B. Loading of Standard Assay tray (maximum of 1)

Secure Assay tray with supplied elastic retention cord


C. Loading of 50ml tubes (maximum of 2)

Note: Two retention cords should be used to secure 50 mL tubes


50 mL tube caps nest securely in debossed area of head


D. Loading of 15 mL tubes (maximum of 2)


E. Loading of 5 mL tubes
(maximum of 2)


F. Loading of 1.5/2 mL tubes
(maximum of 4)

G. Loading of .5 mL tubes
(maximum of 6)

H. Loading of .2 mL tubes
(2 tube strips - maximum of 24)


Vortexer head removal


In the event of a cleaning being needed for the Vortexer, please follow the sequence of removal for the mixing head.

Turn the power **OFF** by placing the 3-position switch into the **OFF** (center) position.


Grasp the installed mixing head attachment with two hands and pull straight up until it releases from the shaft. A technique of lifting the attachment with your finger tips while using your knuckles as a fulcrum will give you the leverage to remove the installed attachment component smoothly. Once Vortexer is cleaned, reinstall the mixing head per instructions.

Pull straight up

Vortexer elastic cord positioning

Care and Handling

Your Vortexer should be given the care normally required for any electrical appliance. Avoid wetting or unnecessary exposure to fumes. The finish can be washed with water (after unplugging) and mild soap or detergents, using a soft cloth or sponge. Avoid exposing any electrical components of the Vortexer instrument to water, cleaning agents, or other liquid! Keep the unit clean by immediately blotting any spills. Replacement parts are available through most laboratory equipment distributors or directly from Heathrow Scientific, LLC.

CAUTION! Unplug from power before cleaning. Do not immerse. Avoid getting any liquid into the Vortexer housing.

Disposal Responsibilities:

The equipment you purchased may contain hazardous substances that could impact the environment. Per regulations on electronic devices in the European Community, you must use the appropriate disposal systems to avoid exposure of these substances to the environment. The disposal systems will reuse or recycle hazardous materials from your equipment responsibly.

The crossed-out wheeled bin symbol invites you to use those systems.

If you need more information on the collection, reuse and recycling systems, please contact your local or regional waste administration. You can also contact us for more information on the environmental performance of our products.

Product Return Information

Important:

Transporting hazardous materials without a permit is a violation of federal law.

Heathrow Scientific® LLC will not accept any product return that is not appropriately cleaned and decontaminated. In the unlikely event of damage to the unit necessitates return, contact Heathrow Scientific® LLC and obtain return authorization before sending your product for service.

Warranty

IN NO EVENT WILL HEATHROW SCIENTIFIC® LLC'S OBLIGATION UNDER THIS WARRANTY EXCEED THE PRICE OF THE PRODUCT.

Limited Warranty

Heathrow Scientific® LLC warrants that Vortexer will be free from defects in workmanship and material for 5 years from the date of purchase.

If you believe that there is a defect in the product, you must, during the warranty period, notify Heathrow Scientific® LLC, provide proof of purchase, and return the product to Heathrow Scientific® LLC with a Return Authorization form. To obtain a Return Authorization form, please call 1-847-816-5070.

If Heathrow Scientific® LLC is properly notified and, after inspection, confirms that there is a defect and the warranty period has not expired, Heathrow Scientific® LLC will repair, modify, or replace the product, at its sole option, at no charge.

OTHER THAN THIS LIMITED WARRANTY, HEATHROW SCIENTIFIC® LLC MAKES NO WARRANTIES, EXPRESS OR IMPLIED, WITH RESPECT TO THE QUALITY OR PERFORMANCE OF THE PRODUCT, INCLUDING BUT NOT LIMITED TO THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR PARTICULAR PURPOSE WHICH ARE HEREBY DISCLAIMED AND EXCLUDED. HEATHROW SCIENTIFIC® LLC WILL IN NO EVENT BE LIABLE FOR ANY LOSS OF USE, LOSS OF PROFITS, CONSEQUENTIAL, SPECIAL, EXEMPLARY OR PUNITIVE DAMAGES.

THIS WARRANTY DOES NOT COVER:

ANY DEFECT OR DAMAGE CAUSED BY IMPROPER OR UNREASONABLE USE OF THE PRODUCT. (HEATHROW SCIENTIFIC® LLC VORTEXER IS DESIGNED FOR USE ONLY BY TRAINED LABORATORY TECHNICIANS. USE BY ANYONE ELSE WILL VOID THIS WARRANTY.)

ANY PRODUCT THAT HAS BEEN, IN HEATHROW SCIENTIFIC® LLC'S SOLE JUDGEMENT, TAMPERED WITH, ALTERED, OR REPAIRED BY ANYONE OTHER THAN HEATHROW SCIENTIFIC® LLC.

ANY PRODUCT THAT IS INOPERATIVE BECAUSE OF: (a) WEAR OCCASIONED BY USE, (b) NEGLIGENCE, (c) ACCIDENT, (d) INCORRECT MAINTENANCE, OR (e) USE UNDER ABNORMAL CONDITIONS OF TEMPERATURE, DIRT OR CORROSION, OR USE WITH ABRASIVE OR CORROSIVE MATERIALS.

ACCESSORY PARTS, SUCH AS RUBBER AND PLASTIC PARTS THAT ARE DAMAGED BY LIQUIDS OR MISUSE.

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All rights reserved. Printed in China.
DIE HAFTUNG VON HEATHROW SCIENTIFIC® LLC ÜBERSTEIGT IN KEINEM FALL DEN KAUFPREIS DES PRODUKTS.
2022 © Heathrow Scientific® LLC